

Supplemental Table 2. Overview of fluorochrome-conjugated antibodies used for phenotypic characterization via flow cytometry (surface staining).

<i>Fluorochrome</i>	<i>Antigen</i>	<i>Clone</i>	<i>Supplier</i>
BUV737	CD45RA	HI100	BD
BUV395	CD4	RPA-T4	BD
BV785	HLA-DR	L243	BioLegend
BV711	CD27	M-T271	BioLegend
BV650	CD279 (PD-1)	EH12.2H7	BioLegend
BV605	TIGIT	A15153G	BioLegend
BV510	CD8	SK1	BioLegend
BV421	CD28	CD28.2	BioLegend
FITC	Vδ2	IMMU389	Beckman Coulter
PE-Cy7	CD39	A1	BioLegend
PE-Texas-Red	CD38	HIT2	BioLegend
PE	Pan γδ	11F2	BD
APC-Cy7	CD19	HIB19	BioLegend
APC-Cy7	CD14	63D3	BioLegend
APC	CD73	AD2	BioLegend
Alexa Fluor 700	CD3	SK7	BioLegend