Supplementary Materials

- This PDF file includes:
- 1. Supplementary Materials and Methods
- 2. Supplementary Figures
- 3. Supplementary Tables

1. Supplementary Materials and Methods

IRI and Cisplatin induced mouse models

Male C57BL/6J mice (8-week-old) were purchased from Chongqing Tengxin Bioscience (Chongqing, China). The construction of IRI-induced AKI mouse model was previously described (Hu et al., 2017). Brifely, mice were placed under general anesthesia for laparotomy, and their bilateral renal infarcts were clipped for 30 minutes by microaneurysm clips, then the mice were fed for another 24 hours for the construction of IRI-induced AKI model. To construct AKI-to-CKD model, bilateral renal ischemia of 35 minutes, followed by reperfusion for 2 weeks was used. To construct cisplatin-induced AKI model, 30 mg/kg of cisplatin was intraperitoneally injected into mice for 3 days as described in previous studies (Li et al., 2018).

To estimate the therapeutic effect of OA, mice were injected with control or OA (20 mg/kg) by tail vein after IRI surgery (Hu et al., 2017) or intraperitoneal cisplatin injection (Li et al., 2018). As reported, the kidney tissues and blood samples were collected (Hu et al., 2017, Li et al., 2018). Serum creatinine (Scr) and blood urea nitrogen (BUN) were determined by using assay kits (Nanjing Jiancheng Bioengineering Institute, China).

Cell Culture

Human proximal tubular epithelial cells (HK-2) were bought from the American Tissue Culture Collection (ATCC, Manassas, VA, USA) and cultured in F12 medium containing 10% fetal bovine serum (Mediatech, Manassas, VA, USA).

RNA Extraction and quantitative PCR (qPCR)

Total RNA was extracted from the mouse kidneys and HK-2 cells according to the manufacturers' instruction with an RNA extraction kit (Beyotime Biotechnology, Shanghai, China), which was reversely transcribed into cDNA with RT Master Mix for qPCR kit (MedChemExpress, Monmouth Junction, NJ, USA). qPCR was carried out by using a SYBR Green qPCR kit (MedChemExpress). The specific primers were listed in Supplementary Table 1 and 2.

Mitochondrial DNA (mtDNA) Copy Number Detection

Using the DNA of HK-2 cells as a template, the mtDNA copy number was detected by qPCR. The target gene mitochondria-encoded cytochrome c oxidase 2 (MT-CO2) and the internal reference glycerol-3-phosphate dehydrogenase (GAPDH) were detected respectively. The specific primers were listed in Supplementary Table 2.

Western Blot

Protein extraction and western blot analysis were performed as previously described (Huang et al., 2020, Huang et al., 2020). The signals were detected by an enhanced

chemiluminescence (ProteinSimple, Los Angeles, CA, USA). The primary antibodies were listed as follows: LC3 (12741; CST), P62 (5114; CST), Fibronectin (26836; CST), and α -SMA (48938; CST), Bcl-2 (3498; CST), HIF1 α (36169; CST), SATT3 (9139; CST), Bax (199677; Abcam), PPAR α (227074; Abcam), BNIP3 (109362; Abcam), COX IV (14744; Abcam).

Cell Viability Assay

After OA treatment of HK-2 cells for 48 hours, Cell Counting Kit-8 (CCK-8) solution (Beyotime Biotechnology, C0037) was added, co-incubating at 37 °C for another 1 hour, and the cell viability was analyzed by using Microplate reader (Thermo Fisher Scientific, Waltham, MA, USA).

Apoptosis Analysis

After the treatment for indicated times, HK-2 cell solutions were harvested for staining with Annexin V-FITC/PI for 15 minutes using an Apoptosis Assay Kit (Beijing Labgic Technology, China), and apoptosis was measured by flow cytometry (Gallios, BD Biosciences).

Reactive oxygen species (ROS) Detection

Adherently cultured cells were incubated with 10 μ M DCFH-DA (Beyotime Biotechnology, S0033S) at 37°C for 20 minutes. Cells were collected for numerical analysis by flow cytometer (Gallios, BD Biosciences).

JC-1 Determination

HK-2 cells were treated with the enhanced mitochondrial membrane potential assay Kit (JC-1) (Beyotime Biotechnology, C2003S) at 37°C for 20 minutes and collected for flow cytometry analysis (Gallios, BD Biosciences).

ATP Measurement

The level of ATP was detected according to the manufacturer's instructions with an enhanced ATP assay kit (Beyotime Biotechnology, S0027).

Transfection of BNIP3 siRNA

To verify whether OA maintained mitochondrial homeostasis by inducing BNIP3 expression, the constructed BNIP3 siRNAs (Biomics Biotechnology, Chongqing, China) were transfected into HK-2 cells by Lipofectamine 3000 (Invitrogen, Carlsbad, CA, USA). The BNIP3 siRNA and Negative control siRNA were listed in Supplementary Table 3.

HBAD-mRFP-EGFP-LC3 assay

HK-2 cells were infected with HBAD-mRFP-EGFP-LC3 adenovirus (Hanbio

Biotechnology, Shanghai, China), according to the manufacturer's protocol. Then, the cells under normoxia or HR conditions were treated with OA for 24 hours, and examined with a confocal microscope (Zesis).

Chromatin immunoprecipitation (ChIP) assay

The procedure for the CHIP experiment was described in our previous study (Huang et al., 2020). Briefly, the sonicated DNA was immunoprecipitated with 2 μ g PPAR α antibody, and then amplified by PCR and qPCR, whose primers were listed in Supplementary Table 4.

Statistical Analysis

All data were presented as mean \pm SEM. Unpaired *t* test or one-way analysis of variance (ANOVA) with Tukey's test were used for statistical analysis. Statistically significant was defined as *P* < 0.05.

REFERENCES


Hu, M. C., Shi, M., Gillings, N., Flores, B., Takahashi, M., Kuro, O. M., et al. (2017). Recombinant alpha-Klotho may be prophylactic and therapeutic for acute to chronic kidney disease progression and uremic cardiomyopathy. *Kidney Int.* 5, 1104-1114. dio: 10.1016/j.kint.2016.10.034

Huang, Y., Wang, S., Zhou, J., Liu, Y., Du, C., Yang, K., et al. (2020). IRF1-mediated downregulation of PGC1alpha contributes to cardiorenal syndrome type 4. *Nat Commun.* 1, 4664. dio: 10.1038/s41467-020-18519-0


Huang, Y., Zhou, J., Wang, S., Xiong, J., Chen, Y., Liu, Y., et al. (2020). Indoxyl sulfate induces intestinal barrier injury through IRF1-DRP1 axis-mediated mitophagy impairment. *Theranostics*. 16, 7384-4400. dio: 10.7150/thno.45455

Li, F., Liu, Z., Tang, C., Cai, J., and Dong, Z. (2018). FGF21 is induced in cisplatin nephrotoxicity to protect against kidney tubular cell injury. *Faseb j.* 6, 3423-3433. dio: 10.1096/fj.201701316R


2. Supplementary Figures



Supplementary Figure 1. OA alleviates cisplatin-induced AKI *in vivo*. (A) Representative micrographs of HE staining of kidney sections from control and cisplatin-induced AKI mice treated with control or OA for 3 days. Scale bars, 1.25 mm (top) and 50 μ m (bottom). (B, C) Effects of OA on serum Scr (B) and BUN (C) (n = 8). (D, E) The mRNA levels of Kim-1 (D) and Ngal (E) were analyzed by qPCR (n = 8). Data are presented as means ± SEM. ****P* < 0.001.



Supplementary Figure 2. OA attenuates mitochondrial injury via restoring mitochondrial energy metabolism *in vivo*. (**A**, **B**) qPCR analysis of the expressions of oxidative phosphorylation and fatty acid oxidation-related genes in kidney tissues from sham and AKI mice treated with control or OA (n = 8). (**C**) Representative images of oil red O staining in kidney tissues from sham and AKI mice treated with control or OA. Scale bar, 50 µm. Data are presented as means \pm SEM. ns: no significance. *P < 0.05, **P < 0.01, ***P < 0.001.



Supplementary Figure 3. OA cannot regulate mitophagy in renal tubular epithelial cells. (A) Western blot analysis of the expressions of LC3, P62 and COX IV in HK-2 cells treated with control or various doses of OA for 24 hours (n = 3). (B) The protein expression of LC3, P62 and COX IV in HK-2 cells exposed to HR followed by incubation with control or OA (10 μ M) were analyzed by Western blot (n = 3). (C) HK-2 cells were exposed to HR after the infection with HBAD-mRFP-EGFP-LC3 adenovirus and treated with control or OA (10 μ M) for 24 hours for confocal microscopy observation. Bar = 10 μ m. (D) The protein expressions of LC3, P62 and COX IV in the kidney tissues from AKI mice treated with control or OA (20 mg/kg) were analyzed by Western blot (n = 3). Data are presented as means ± SEM. ns: no significance. ***P* < 0.01, ****P* < 0.001.

3. Supplementary Tables

Gene(mouse)	Primers	Product length	
β -actin	Forward: 5'-TGTTACCAACTGGGACGACA-3'	1 <i>65</i> h.:	
	Reverse: 5'-GGGGTGTTGAAGGTCTCAAA-3'	165 bp	
CROT	Forward: 5'-GGCAAAAAGCTCACCAAGGA-3'	63 bp	
	Reverse: 5'-CTGAAGAGCGAGCTGAATAAAGG-3'		
Cpt1b	Forward: 5'-TCTAGGCAATGCCGTTCAC-3'	00.1	
	Reverse: 5'-GAGCACATGGGCACCATAC-3'	99 bp	
ACADM	Forward: 5'-AGGATGACGGAGCAGCCAATGA-3'	120 h-	
	Reverse: 5'-GCCGTTGATAACATACTCGTCAC-3'	120 bp	
ACSL-1	Forward: 5'-AAATTCTCGCGGTGTCTCCA-3'	1501	
	Reverse: 5'-CACTCGGGGGTCTACGTTTCA-3'	150 bp	
Atp5a1	Forward: 5'-GCCCTCGGTAATGCTATTGA-3'	200 1	
	Reverse: 5'-GCAATCGATGTTTTCCCAGT-3'	209 bp	
Kim-1	Forward: 5'-ACATATCGTGGAATCACAACGAC-3'	114 bp	
	Reverse: 5'-ACTGCTCTTCTGATAGGTGACA-3'		
N	Forward: 5'-GCAGGTGGTACGTTGTGGG-3'	95 bp	
Ngal	Reverse: 5'-CTCTTGTAGCTCATAGATGGTGC-3'		
Cytc	Forward: 5'-ACCAAATCTCCACGGTCTGTT-3'	126 h.,	
	Reverse: 5'-GGATTCTCCAAATACTCCATCAG-3'	126 bp	
N. J. C.	Forward: 5'-ATCCCTTACCCTTTGCCACT-3'	100 h.	
Ndufa	Reverse: 5'-CCGTAGCACCTCAATGGACT-3'	198 bp	
	Forward: 5'-TGTGAGACCGTGCTAATGGA-3'	2171	
Ndufv	Reverse: 5'-CATCTCCCTTCACAAATCGG-3'	217 bp	
	Forward: 5'-TCAGGGTACCACTACGGAGT-3'	2061	
PPARa	Reverse: 5'-CTTGGCATTCTTCCAAAGCG-3'	206 bp	

Supplementary Table 1. The primer sets for mouse

Gene(human)	Primers	Product length	
β-actin	Forward: 5'-GTGAAGGTGACAGCAGTCGGTT-3'	157 bp	
	Reverse: 5'-GAAGTGGGGTGGCTTTTAGGA-3'		
Nd.f.	Forward: 5'- ATGAAGGTGACAGCGTGAGG-3'	100 bp	
Ndufv	Reverse: 5'- TTCTTGGGTGCTGTCCGC-3'	100 bp	
Marka	Forward: 5'-GCGCTGTTGCCACTGTT-3	102 hp	
Ndufa	Reverse: 5'-AGTCGGGCACGTCCTTC-3	103 bp	
A 4m 5 m 1	Forward: 5'-ATTTCTCAGTCTACGCCGCA-3'	147 -	
Atp5a1	Reverse: 5'-AGACACGCCCAGTTTCTTCA-3'	147 bp	
ACADM	Forward: 5'- GGCCGTGACCCGTGTATTAT-3'	77 ha	
ACADM	Reverse: 5'-CTGCAGCATCGCCCGAA-3'	77 bp	
Critlh	Forward: 5'- AGGATCTGGGCTTACCTAGAG-3'	144 hm	
Cpt1b	Reverse: 5'-ACTTGCCCACCATGACTTGA-3'	144 bp	
CROT	Forward: 5'- GGACACCCTGGTTGTTGCTA-3'	155 hr	
CKOI	Reverse: 5'-AACATCTTTTGCTGCCGCTC-3'	155 bp	
Criste	Forward: 5'-CGTTGTGCCAGCGACTAAAAA-3'	100 ha	
Cycs	Reverse: 5'-TTGCCTCCCTTTTCAACGGT-3'	109 bp	
	Forward: 5'-GCTCCCAGACACCACAAGAT-3'		
BNIP3	Reverse: 5'-TGAGAGTAGCTGTGCGCTTC-3'	225 bp	
	Forward: 5'-TAAAGACAGATGGGAGGAGACCC-3'	170 h -	
ACSL-1	Reverse: 5'-ACGTACTGTCGGAAGTCAACC-3'	170 bp	
0040	Forward:5'-GCTATCATTACGGAGTCCACG-3'	001	
PPARa	Reverse: 5'-TCGCACTTGTCATACACCAG-3'	88 bp	
MT-CO2	Forward: 5'-CAAACCTACGCCAAAATCCA-3'	1 < 41	
M1-CO2	Reverse: 5'-GAAATGAATGAGCCTACAGA-3'	164bp	
CADDII	Forward: 5'-TGACAACAGCCTCAAGAT-3'	1051-	
GAPDH	Reverse: 5'-GAGTCCTTCC ACGATACC-3'	195bp	
	Forward: 5'-CATAAAGCAGGCAGAACTCA-3'	1071	
TFAM	Reverse: 5'-GCAACTAAAACCAGCAATCAC-3'	107bp	

Supplementary Table 2. The primer sets for human

Supplementary Table 3. The sequence sets for siBNIP3 and siNC

Gene(human)	Sequence		
siBNIP3	Sense:	5'- GGAACACGAGCGUCAUGAA -3'	
SIDIVII J	Antisens	e: 5'- UUCAUGACGCUCGUGUUCC -3'	
	Sence:	5'- UUCUCCGAACGUGUCACGUTT -3'	
siNC	Antisens	e: 5'- ACGUGACACGUUCGGAGAATT -3'	

Gene (human)	Primers	Product length	
BNIP3	Forward: 5'- GGCAAATACCTCTTGATGTGCATAAAC -3'	146 bp	
	Reverse: 5'- ACCTATTCGAGGGAGTCAACTAAAAGT -3'		

Supplementary Table 4. The primer sets for ChIP