

## **Supplemental Figure 2:**

Protein phosphorylation by native or aggregated antibodies. (A) Syk phosphorylation is concentration-dependent in moDC treated by infliximab HSL2 aggregates. (B) Syk, ERK1/2 and Akt phosphorylation in moDC treated with rituximab or adalimumab and aggregates. Phosphorylation was detected following 15 and 30 min of stimulation. Results are expressed as fold change over PBS control. Results from 12 independent experiments are shown. \*p<0.05; \*\*p<0.01; \*\*\*p<0.001.